

# HOBOWare

From labwiki

## Uploading data

1. Open HOBOWare Pro Trial
2. Connect the HOBOWare via its USB port. When the device is successfully connected '1 device connected' will appear in the bottom right hand corner of the window.
3. Click on the 'Readout device' icon, found in the top left hand corner (second icon).
4. Change the temperature settings from Fahrenheit to Celsius.
5. Click 'plot'.
6. Save the data.
7. Export data as an excel file by clicking the 'Export Table Data' icon on the tool bar. Save the csv. file generated.

Retrieved from "<http://ecothery.uwaterloo.ca/labwiki/index.php?title=HOBOWare&oldid=22>"

- 
- This page was last modified on 30 April 2014, at 10:55.
  - This page has been accessed 33 times.